
Bedienungsanleitung
Instruction manual
Manual de instrucciones

Pyro Light Optima **Pyro Light Optima XL**

DE

EN

ES

Enjoy mobility.

Inhalt

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Contents

1.	Introduction	56
1.1.	Labels and symbols	56
2.	Important safety information	57
2.1.	General safety information	57
2.2.	Risk of tipping over	60
2.3.	Safety in vehicles	60
2.4.	Use on roads	60
2.5.	Brake	61
2.6.	Intended use	61
2.7.	Indications	61
2.8.	Contraindications	62
2.9.	Declaration of conformity	62
2.10.	Liability	62
2.11.	Service life	63
3.	Product and scope of delivery	63
3.1.	Checking the delivered goods	63
3.2.	Scope of delivery	63
3.3.	Product's identification label and serial number	64
3.4.	Overview	65
4.	Assembly/Adjustments	66
4.1.	Folding and unfolding	66
4.2.	Legrests	67
4.2.1.	Swing-away legrests	67
4.2.2.	Adjusting the lower leg length	68
4.3.	Angle-adjustable footplates	68
4.4.	Elevating legrests (optional)	68
4.5.	Amputee legrests (optional)	69
4.6.	Sideguards	69
4.6.1.	Standard sideguard	69
4.7.	Seat	70
4.7.1.	Seat depth	70
4.7.2.	Seat height	71
4.7.3.	Seat angle	72
4.8.	Wheelbase	73
4.9.	Wheels and tyres	74

4.9.1.	Quick-release axle	74
4.10.	Brakes	74
4.10.1.	Operating the parking brake	75
4.10.2.	Adjusting the parking brakes	76
4.11.	Operating the parking brake with the brake lever extension	76
4.12.	Drum brake (optional)	77
4.13.	Backrest	78
4.13.1.	Back sling	78
4.13.2.	Backrest height	78
4.14.	Push handles	79
4.15.	Tipping aids	79
4.16.	Anti-tip wheels (optional)	79
4.17.	Reflectors	80
5.	Accessories	80
5.1.	Lap belt (optional)	80
5.2.	Therapy table (optional)	80
5.3.	Togo (optional)	80
5.4.	Crutch holder (optional)	81
5.5.	Spokeguard (optional)	81
5.6.	Headrest	81
5.7.	One-hand operation (optional)	81
5.8.	One-hand brake (optional)	81
5.9.	Handrim cover (optional)	82
5.10.	Infusion & oxygen bottle holder (optional)	82
5.11.	Fixing a bag (optional)	82
6.	Use	83
6.1.	Side transfer	83
6.2.	Frontal transfer	84
6.3.	Going up or downstairs or a high obstacle	84
6.4.	Driving the wheelchair	85
7.	Technical data	86

8.	Transport	88
8.1.	Transporting the wheelchair	88
8.2.	Transportation in vehicles	88
9.	Maintenance instructions for the user	89
9.1.	Cleaning	89
9.2.	Disinfection	90
9.3.	Checks before driving	91
9.3.1.	Checking the brakes	91
9.3.2.	Checking the tyres	91
10.	Repairs	91
11.	Maintenance information	92
12.	Passing the wheelchair on to another user	92
13.	Storage / shipping	92
14.	Disposal	93
15.	Warranty	94

1. Introduction

Dear user,

You have chosen a high-quality wheelchair from Bischoff & Bischoff. Thank you for trusting in us.

This user manual contains all the information necessary to operate the Pyro Light Optima / Pyro Light Optima XL wheelchairs.

Your wheelchair was designed for both indoor and outdoor use.

Before using your new wheelchair for the first time, please read and follow the information in this user manual carefully. This is an integral and necessary part of the wheelchair. Keep the user manual within easy reach and pass it on when transferring the wheelchair.

For users who are visually impaired, this document is available as a pdf file on our website www.bischoff-bischoff.com

Note!

Despite thorough research and editing, incorrect information may occur in our user manual. We reserve the right to make technical modifications. Images may show accessories. You can find the current version of the user manuals in the download area of our website:

www.bischoff-bischoff.com

1.1. Labels and symbols

Important! Indicates particularly useful information in the relevant context.

Attention! Means particular safety-related Notes. Observe the user manual!

2. Important safety information

2.1. General safety information

- To avoid falls and dangerous situations, you should first practise using your new wheelchair on level, manageable terrain or in an obstacle-free area. An accompanying person is recommended in this case.
- Only drive on level, solid ground.
- Be aware of the danger of jamming between moving parts when adjusting, using and maintaining the wheelchair.
- The accompanying person should be physically and mentally capable of driving a wheelchair.
- The driver and attendant must not be under the influence of alcohol or drugs.
- Before using the wheelchair, please check that all attached parts are properly fastened.
- Only lift the wheelchair by firmly attached parts (the push handles or the frame parts).
- Make sure that the quick-release axles are correctly seated in the axle sockets, otherwise there is a risk of injury.
- Keep clothing and body parts away from gaps between moving parts.
- Please note that the risk of tipping over is increased if there is a shift in balance (e.g. due to strong upper body movements or when overcoming obstacles).
- When negotiating obstacles (steps etc.), an attendant must be called in. Use the tipping aid (see chapter 4.15) to tilt the wheelchair accordingly.
- Stairs and larger obstacles may only be overcome with the help of at least two accompanying persons. Use the lifting points for this (see chap. 6.3). If facilities such as ramps, lifts or stair-climbing aids are available, these must be used (Fig. A).
- When driving on a slope, please note that, depending on the axle position, there is a risk of tipping over from just a few degrees.

Pyro Light Optima | Lightweight wheelchair Pyro Light Optima XL

- Avoid collisions of any kind! Avoid driving into an obstacle (step, kerb) without braking or jumping down steps (Fig. B).
- Replace the tyres as soon as the tread is worn down.
- When getting in or out of the wheelchair, apply both parking brakes.
- Only use the parking brake for parking, it must not be used to slow down the ride.
- Avoid driving on uneven, loose surfaces.
- In road traffic, the road traffic regulations must be observed (see chap. 2.4).
- The seat and upholstery materials meet the requirements of resistance to flammability according to EN 1021- 2. Nevertheless, do not expose upholstery or other parts to sources of fire - such as cigarettes.
- The wheelchair may only be used to transport one person. The wheelchair must not be used to transport loads or more than one person.
- When braking from fast travel or on long downhill stretches, using the handrims, fingers and palms heat up. Caution, danger of burns! (Fig. C).
- Avoid fast travel and downhill stretches.
- Note that upholstered parts exposed to direct sunlight heat up and can cause skin injuries if touched. Therefore, cover these parts or protect the wheelchair from sunlight.

Protect the wheelchair from sunlight to avoid burns caused by heated components.

- Also note that the frame and upholstered parts can cool down considerably in winter. If possible, do not park the wheelchair outside in cold weather.
- Observe the storage and operating conditions.

Lightweight wheelchair | Pyro Light Optima Pyro Light Optima XL

- When transferring between the wheelchair and the bed, always observe the following:
- apply both parking brakes.
- swing the side part of the entry/exit side downwards.
- swing the foot plates (when getting in/out) up or away.
- Reporting incidents.

If you, as the operator, user or their relatives, detect possible defects or functional limitations, please contact an authorised dealer immediately. Operators, users or their relatives should inform the dealer who supplied the dealer who provided you with the product of any suspected serious incidents which directly or indirectly have had or could have one of the following and by which they are affected. You can also report this to the competent higher federal authority report

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of the health of the state of health of a patient, user or other persons,
- a serious risk to public health.

2.2. Risk of tipping over

If fitted, the anti-tip wheels which are fitted on both sides (see Chapter 4.16) will largely prevent the wheelchair from tipping over backwards. However please note that driving with anti-tip wheels is only possible in a limited way on ground with different levels, lift platforms and ramps. Make sure there is sufficient gap at the top, bottom and sides.


2.3. Safety in vehicles

To see whether your wheelchair is permitted for use as a seat in a Vehicle for Transporting Individuals with Impaired Mobility or not, you can refer to the sticker on the identification label on the wheelchair.

- Wheelchair is designed for use as a seat in a Vehicle for Transporting Individuals with Impaired Mobility (Fig. D).
- Wheelchair is not designed for use as a seat in a Vehicle for Transporting Individuals with Impaired Mobility (Fig. E)

If your wheelchair is not designed to be used as a seat in a vehicle, under no circumstances must it be used as such. It does not fulfil the requirements in compliance with ISO 7176-19. Violation can cause the most serious injuries resulting in death in the case of an accident (Fig. B).

2.4. Use on roads

The wheelchair is designed for indoor and outdoor use. Please note that when you drive on public roads you must comply with road traffic regulations. Do not endanger other users through reckless driving, this applies particularly on pavements.

2.5. Brake

Braking during use:

When in motion, brake the wheelchair using the handrims. Please note that the hands may become hot when doing this (see Chapter 4.10).

Parking brake/parking:

As soon as the wheelchair has stopped, push both brake levers forward (Fig. 23). The wheelchair is securely braked when it can no longer be pushed and both brake levers are pushed forward to the stop position. Pull the lever towards you to release the brake (see Chapter 4.10.1).

Drum brake:

As an option, the wheelchair can be fitted with a drum brake which can only be operated by an attendant using the brake lever on the push handles.

2.6. Intended use

The wheelchair is designed to increase the mobility only for adult persons and those who have difficulty walking with the indications specified (see Section 2.7). The maximum payload of

- 125 kg for the Pyro Light Optima and
 - 170 kg for the Pyro Light Optima XL
- must not be exceeded.

Please note:

We can only provide a warranty if the product is used under the conditions specified and for the intended purposes.

2.7. Indications

Provision of a wheelchair is intended for adult persons who are unable to walk or who have severely impaired mobility, e.g. through:

- paralysis
- loss of limbs
- limb defects/deformity

- joint contractions (not on both arms)
- joint damage (not on both arms)
- other illnesses

2.8. Contraindications

Provision of wheelchairs is not suitable for adult persons who have:

- no capability to sit
- joint damage/contractions on both arms
- perceptual disorders
- severe equilibrium disorders
- loss of limbs on both arms
- reduced and insufficient vision
- severe limitations of cognitive abilities

2.9. Declaration of conformity

We, as Bischoff & Bischoff GmbH, declare under sole responsibility, that the Pyro Light Optima wheelchair meets all requirements of MDR 2017/745, which are applicable.

2.10. Liability

We can only provide a warranty if

- the product is used under the conditions specified and for the intended purpose.
- Modifications, enhancements, repairs and maintenance work is carried out only by those individuals who have been authorised by us to do so, and
- the wheelchair is used in accordance with the instructions for use.

2.11. Service life

The expected product service life is up to five years. Prerequisites for this are proper use as well as compliance with the safety, care and maintenance information in accordance with this user manual.

3. Product and scope of delivery

3.1. Checking the delivered goods

All products from Bischoff & Bischoff GmbH undergo an appropriate final inspection at our factory and the CE mark is affixed (see Chapter 3.3).

The wheelchair will be delivered in a special cardboard box. After unpacking, retain the cardboard box if possible. It can be used for possible later storage or return of the product.

Delivery and induction is usually carried out by a qualified medical equipment retailer. When sending by train or shipping company, the goods must be checked immediately for transport damage in the presence of the forwarding agent.

Check that the delivery is complete and undamaged (see Chapter 3.4). In the case of anomalies or damage, contact our Customer service. You can find the contact details on the back cover of this user manual.

3.2. Scope of delivery


On receiving the goods, please immediately check that the delivery is complete. The delivery consists of:

- packing,
- wheelchair (preassembled)
- legrests
- this user manual
- any necessary accessories

3.3. Product's identification label and serial number

The product's identification label (Fig. 1) and serial number are located on the cross brace.

- A** Manufacturer's logo
- B** Manufacturer information
- C** MD Medical product
- D** CE mark
- E** Attention! - read the User manual
- F** Important! - read the User manual
- G** Model name
- H** Model number
- I** Max. payload
- J** Date of manufacture
- K** UDI number
- L** Serial number


Depending on the model number, the information may vary. Refer to your product's identification label for the information which relates to your product.

3.4. Overview

(see Fig. 2)

- A** height-adjustable push handles
- B** Adjustable backrest and upholstery
- C** Castor
- D** Sideguard
- E** Seat sling
- F** Drive wheel
- G** Quick-release axle
- H** Parking brake
- I** Handrims
- J** Legrest with heel strap


4. Assembly/Adjustments

Your new Bischoff & Bischoff wheelchair will be delivered in a Bischoff & Bischoff original cardboard box, assembled and folded.

In order to avoid transport damage, detachable attachments are delivered separately with the product.


Refer to the technical data to see which options are possible on your wheelchair (see Chapter 5).

4.1. Folding and unfolding

Unfolding

- Position the folding wheelchair next to you (Fig. 3).
- Tip the wheelchair towards you until the weight is off one of the drive wheels.
- Push the seat tubes apart and push them down with the flat of your hand, until they lock into place in the bracket on the frame (Fig. 4).
- The seat pad must be completely unfolded.
- Apply the left and right parking brakes on the drive wheels.
- Fit the legrests.
- You can now sit down.

Folding

- Position the folding wheelchair next to you.
- Flip the footplates up (Fig. 5).
- Remove the seat cushion if one is fitted, take hold of the front and rear of the seat upholstery and pull it up.

When folding and unfolding, watch your fingers as there is a risk of pinching.


4.2. Legrests

4.2.1. Swing-away legrests

You can flip the footplates up separately to the side and adjust to the lower leg length of the user. You can remove the legrests, swing them away and refit them without needing to use tools.

Attach legrests

- Hold the top of the legrest so that it is pointing outwards.
- Guide the legrest into the hole in the frame using the guide bolts (Fig. 6).
- Swing the legrest to the front until you hear it clearly click into place (Fig. 6).
- Flip the footplate down if necessary.

Swing the legrests away and remove them



In order to make it easier to transfer to and from your folding wheelchair, you can swing the legrests inwards or outwards, without having to remove them.

- Flip the footplate up if necessary.
- Pull the handle upwards and swing the legrest outwards (Fig. 7/8).
- In the outward position, the legrest can be pulled out by pulling it upwards (Fig. 8).


Before setting off, fit the legrests as described above.

Never stand on the footplate and when swinging it away, watch your fingers as there is a risk of pinching.

4.2.2. Adjusting the lower leg length


Loosen the screws, then adjust the legrest to the required length and re-tighten the screws (Fig. 9).

4.3. Angle-adjustable footplates

The angle of the footplates can be adjusted:

- Unscrew the hexagon screw at the sprocket (see Fig. 10). Move the footplate to the required position by rotating the sprocket.
- Re-tighten the sprocket screw (see Fig. 10).

4.4. Elevating legrests (optional)

Assemble, disassemble and swing them away as for the standard legrests.

To adjust, undo the clamp lever, then adjust the legrest to the required length and then tighten the clamp lever

The horizontal adjustment works as follows:

Upwards:

- By lifting the footplate receiver tube you can get the legrest into the required position.

Downwards:

- Pull locking lever upwards.
- When doing this, hold the legrest firmly and move to the required position.
- Release locking lever; legrest locks into place automatically.

Never stand on the footplate and when swinging it away, watch your fingers as there is a risk of pinching.

For transfer, the legrest must be swung away or removed. Never sit on the legrest. Risk of tipping over!

4.5. Amputee legrests (optional)

Assemble, disassemble and swing them away as for the standard legrests, adjustment works as for the horizontally-adjustable legrests (see chapter 4.4).


For transfer, the legrest must be swung away or removed. Never sit on the legrest. Risk of tipping over!

4.6. Sideguards

4.6.1. Standard sideguard

- In order to be able to transfer more easily, both sideguards with armrests can be swung backwards and removed.
- Pull the locking lever gently and at the same time, swing the sideguard to the back (Fig. 11).
- Once the front side of the sideguard is released and the fixing bolt at its rear has been pulled out, the sideguard can be removed by an upwards movement.
- To lock the sideguard, flip it forwards and click it into place on the armrest using gentle pressure.


When swinging the sideguards away, make sure that clothing or other items cannot be trapped. Before using for the first time, check that the sideguards are firmly in place!


Never use the sideguards to carry the wheelchair!

Height adjustment

- Pull the locking lever to move the armrest to the required position (Fig. 12).
- When you have reached the position required, then release the lever. The lock clicks into place.
- The height has been correctly selected if the upper arm and forearm are approximately at right angles to one another when the forearm is resting on the armrest.


Depth adjustment

- Press the knob on the side of the armrest and push the armrest to the required position (Fig. 13).
- Release the knob and move the armrest gently forwards or backwards until it clicks into place.

4.7. Seat

4.7.1. Seat depth

Attention: risk of tipping over! If you move the drive wheel forwards or the backrest backwards, then the risk of tipping over increases significantly. Anti-tip wheels or a wheelbase extension are required. Check the settings first with an assistant. Only use settings which you consider to be safe.

be adjusted in a range of 10 centimetres. First remove the bracket for the sideguard and backrest by undoing the fixing screws (Fig. 14). Then remove the hexagon screw at the back of the seat sling (Fig. 15). Move the adapter to the desired position. Then pull the rear part of the seat upholstery back and connect it to the front seat upholstery using the hook and loop fastener (Fig. 16). Then re-tighten the screw of the adapter.

Fix the backrest back onto the wheelchair frame. Position the backrest tube as close to the seat as possible. Then fit the bracket so that the sideguards click into place properly.

The seat depth can also be adjusted by the position of the backrest.

For safety reasons, if the backrest is set to the back, then select the largest wheelbase possible or use a wheelbase extension.




4.7.2. Seat height


You can adjust the seat height to meet your requirements in just a few steps.

Drive wheel height

Remove the drive wheel (see Chapter 4.9.1) to make this task easier. To adjust the height of the drive wheel and thereby raise or lower the entire wheelchair, you must first remove the drive wheel axle sleeve. You can undo the fixing nuts (Fig. 17) of the axle sleeve using a spanner. You can now remove the sleeve and adjust to the required height. Re-tighten the nuts and fit the drive wheel.

When the drive wheel height is changed, the parking brake must be adjusted (Chapter 4.10.2).


Height of the front wheel

The front wheel can be adjusted to three heights. Undo the hexagon lock nut (Fig. 18) and pull out the wheel axle. Mount the front wheel at the required height and re-tighten the lock nut.

In addition the castor bearing bracket can be rotated to achieve even greater adjustment.


Angle of the steering wheel axle

By changing the height of the drive wheel or front wheel, the angle of the steering wheel axle may change. However this should always be at right angles (90°) to the floor (Fig. 19).

First undo the lock by loosening the two Allen screws on the inner side (Fig. 20). You can then adjust the angle using the two setscrews on the under side (Fig. 20).

Then re-tighten the two fixing screws on the inner side.

4.7.3. Seat angle

If you tilt the seat surface of the wheelchair to the back, you will therefore be sitting more stable in the wheelchair. To tilt the seat backwards, you must fit the quick-release axle in a higher position. If you want to tilt the seat forwards, fix the quick-release axle in a lower position.

4.8. Wheelbase


Attention: risk of tipping over! If you move the drive wheel forwards or the backrest backwards, then the risk of tipping over increases significantly. Anti-tip wheels or a wheelbase extension are required. Check the settings first with an assistant. Only use settings which you consider to be safe.

The wheelchair's stability and steadiness is improved by extending the wheelbase.

Conversely, if the wheelbase is reduced, the manoeuvrability is improved. Reducing the wheelbase is only recommended for experienced riders.

To adjust the wheelbase, undo the four fixing screws which fix the wheel support (Fig. 21). The wheelbase can adjusted in four positions. Then refit the wheel support in the recesses provided for this purpose. Make sure that all screws are firmly tightened.

If the wheelbase is changed, then it is essential to adjust the knee lever brake as well!

4.9. Wheels and tyres

Your wheelchair is fitted with PU tyres as standard.

If the drive wheel or castor is worn or defective, ask your dealer to service them!

4.9.1. Quick-release axle

As standard the wheelchair has quick-release axles. This means you can easily remove and refit the drive wheels.

- Push the lock of the quick-release axle using your thumbs and pull the wheel off (Fig. 22).
- When refitting, simply put the wheel onto the receiver. When doing this, push the locking knob in.

The click you will hear from the quick-release axle will indicate that the wheel is securely mounted. After fitting, make sure that the wheels are tight!

4.10. Brakes

As standard, the wheelchair is fitted with parking brakes to be operated by the wheelchair user.

The parking brakes must not be used as operational brakes, as there is a risk of death due to the risk of falling!

Please note that there is a risk of burns if your hands heat up when braking using the handrims!

4.10.1. Operating the parking brake


The parking brake must only be used for parking. To slow down while driving, you must only use the drive wheel handrims or the optional drum brakes for attendant.


- When braking using the handrims, hold the highest point. First let the handrims slide through your hands
- Slowly increase the braking force by exerting gripping force on the handrims, until you come to an abrupt stop.
- As soon as the wheelchair has stopped, push both brake levers forward (Fig. 23). The wheelchair is securely braked when it can no longer be pushed and both brake levers are pushed forward to the stop position.

Pull the lever towards you to release the brake.

Never use the brake lever as a support - e.g. when transferring or standing up -, otherwise there is a risk of tipping over and falling! The lever could break off!


The parking brake - with or without brake lever extension - only works on one drive wheel and must therefore always be used on both sides!

4.10.2. Adjusting the parking brakes

Adjusting the brakes must only be carried out by trained personnel. Please bear in mind that if the brakes are set incorrectly this can have life-threatening consequences.

- Mount the brake in a position which is operable and offers secure braking in all driving situations.
- To adjust the gap, undo the two socket screws on the brake bracket (Fig. 24).
- Then push the brake system to the correct position along the bracket.
- Finally re-tighten the two socket screws.
- Make sure that the brakes click into place when operating the lever and that they don't snap back.


4.11. Operating the parking brake with the brake lever extension

The procedure for operating the parking brake with the brake lever extension is the same as operating the parking brake (see chapter 4.10.2).

Folding the brake lever extension in

- The brake lever can be folded in when driving or transferring to the side.
- Pull the plastic grip out of the lever (Fig. 25).
- Fold the plastic grip down (Fig. 26).

- Before setting off, the brake must be returned to its operational position by following the procedure in reverse order.

Never use the brake lever as a support - e.g. when transferring or standing up -, otherwise there is a risk of tipping over and falling! The lever could break off!

The parking brake - with or without brake lever extension - only works on one drive wheel and must therefore always be used on both sides!


4.12. Drum brake (optional)


The drum brake does not replace the parking brake.

Adjusting the brakes must only be carried out by trained personnel. Please bear in mind that if the brakes are set incorrectly this can have life-threatening consequences.

As an option, the wheelchair can be fitted with a drum brake which can only be operated by an attendant using the brake lever on the push handles (Fig. 27).

Adjusting the drum brake

Basic setup for the drum brake is carried out on the cable clamping screw. Fine adjustment is carried out using the adjusting nut on the brake cable. Unscrewing the adjustment nut will increase the braking effect (Fig. 28).


- Drum brakes only work on one drive wheel at a time!
- Drum brakes are also subject to natural wear. If the braking effect decreases, the wheelchair must no longer be used. In order to avoid this happening, for your own safety you should check that the brakes are working properly each time before setting off!
- Once you have finished the adjustment, re-tighten all screws and carry out a brake test!

4.13. Backrest

In order to ensure a comfortable seating position, the rigidity of the back sling can be adjusted. You can also adjust the backrest height.

4.13.1. Back sling

The more taut the back sling, the more rigid the backrest will feel. First remove the back cushion which is fixed in place with a hook and loop fastener. You can then change the tension of the individual straps so that they are more or less taut. In order to relieve tension on certain areas of your back, you can make the straps there slightly looser. You must now refix the back padded cover with the hook and loop fastener.


4.13.2. Backrest height

In order to be able to adjust the backrest height, we recommend removing the sideguards and back padded cover first. You can then undo the screw fitting on both sides (Fig. 29) and move the frame to

the required height. Then refit the screws in the recesses provided for this purpose and tighten them by hand.


4.14. Push handles

The height of the push handles can be adjusted to make it more comfortable for the attendant. Undo the lever screw (Fig. 30) on the push handle and move to the required height. Then re-tighten the lever screw to lock the push handle in place.


4.15. Tipping aids

In order to make it easier for an attendant to negotiate obstacles such as kerb edges, tipping aids can be used. Hold the wheelchair firmly with both hands on the push handles. Press one foot onto one of the tipping aids and at the same time push the hands down (Fig. 31). Tipping may be restricted by the anti-tip wheels.

4.16. Anti-tip wheels (optional)

The anti-tip wheels prevent the wheelchair from tipping over backwards. There must be a gap between the rollers of the anti-tip wheels and the ground of approx. 30-50 mm (Fig. 32). Anti-tip wheels are particularly recommended for inexperienced wheelchair users or amputees.

For safety reasons, always use the anti-tip wheels as a pair!

4.17. Reflectors

As standard your wheelchair is fitted with reflectors on the back and spokes. These must not be covered.

5. Accessories

5.1. Lap belt (optional)

You can order a lap belt as an option and this will secure users who do not have the required stability in the wheelchair.

- Remove the sideguard or swing it away (see Chapter 4.6).
- On both open ends, wrap the belt once round the sideframe tube.
- Then fix the belt strap with the buckle.
- You can adjust the belt length using the sliding tabs on the belt.
- Use the fastening clip to open and close the belt.

The lap belt must always be fastened and closed when using the wheelchair.

5.2. Therapy table (optional)

The table is simply pushed into the armrest from the front. To remove the therapy table, simply pull it back out.

5.3. Togo (optional)

An electrical Togo pushing & braking aid can be fitted to your wheelchair.

If required, see the Togo pushing & braking aid user manual for all relevant information on assembling and use.

5.4. Crutch holder (optional)

Put the crutch in the crutch holder and then fix it to the backrest using the retaining strap.

Never step on the crutch holder and never use it as a tipping aid.

5.5. Spokeguard (optional)

The spokeguard is fixed to the spokes using the cable ties.

5.6. Headrest

The headrest is screwed onto the push handles. After loosening the clamping lever or star screw, you can adjust the position of the headrest. Then tighten the screws firmly.

5.7. One-hand operation (optional)

With one-hand operation you have both handrims on one wheel, so that you can operate both drive wheels with one hand.

The option can be fitted either left or right.

5.8. One-hand brake (optional)

With the one-hand brake you have a parking brake where you can brake both drive wheels by pressing only one brake lever.

The option can be mounted on the left or right.

5.9. Handrim cover (optional)

Put the handrim cover over the handrims in order to achieve better grip on the handrim, for example if your hand strength is limited.

With better grip, there is a risk of sharp braking.
Attention: risk of falling!

5.10. Infusion & oxygen bottle holder (optional)

The infusion & oxygen bottle holder is pushed onto the frame tube behind the anti-tip wheels and fixed to the backrest of your wheelchair. Attention - with increased load from the oxygen bottles, the risk of your wheelchair tipping over backwards will also increase. We strongly recommend the use of anti-tip wheels.

5.11. Fixing a bag (optional)

If you want to fit a bag to your wheelchair, hang the straps of your bag over the back tubes of your wheelchair.

Please note that the reflectors must not be covered. Bags must not interfere with the spokes of the drive wheels. Attention: risk of falling!

Attention - with increased load in the bag, the risk of your wheelchair tipping over backwards will also increase. We strongly recommend the use of anti-tip wheels.

6. Use

6.1. Side transfer

Never stand on the footplates, there is a Risk of tipping over!

- Make sure that the sides of the wheelchair and the current or future seat surface are as close to each other as possible.
- If you are transferring from another wheelchair or indoor frame, secure these by applying the parking brakes.
- Secure the wheelchair using the parking brake (see Chapter 4.10).
- If necessary, swing the armrest backwards or remove the sideguard (see chapter 4.6).
- Flip both footplates up, swing them to the side and remove them (see chapter 4.2).
- Slide sideways onto the other seat surface. When doing this sit as far back as possible.
- Move the legrests back into their original position.
- If necessary lock the sideguards back in place.

6.2. Frontal transfer

Never stand on the footplates, there is a Risk of tipping over!

- Remove the legrests (see Chapter 4.2).
- Make sure that the front of the wheelchair and the current or future seat are as close to each other as possible and apply the parking brakes on both sides (see Chapter 4.10).
- If you are transferring from another wheelchair or indoor frame, secure these by applying the parking brakes.
- Now slide onto the seat surface by rotating your body.
- Then hang the legrests back into the initial position and flip the footplates down.

6.3. Going up or downstairs or a high obstacle

Only the frame and the push handles must be used as grip points.

The helpers only grip on the firmly fixed frame parts (Fig. 33). In order to go up stairs, the rear helper pulls the wheelchair backwards up the steps by holding onto firmly-mounted handles. The lower helper grips on the front frame tubes and stabilises the wheelchair. When doing this he pushes the drive wheels of the wheelchair onto the steps. When descending, the lower helper brakes the chair by pressing the wheelchair straight into the steps. The upper helper holds the wheelchair by the push handles, keeps it secure and holds it in the right position. The wheelchair should roll down step by step.

Two attendants must help with obstacles and stairs with more than three steps.

6.4. Driving the wheelchair

- Sit in your Bischoff & Bischoff wheelchair. Make sure that you sit as far back as possible.
- Propulsion is achieved using the handrims on the drive wheels. Place your hands on the handrims using your thumbs and bent index fingers.
- From this position, push the wheelchair with both hands (Fig. 34).
- Turning right: hold the right handrim tight and only push with the left one.
- Turning left: push in the opposite way.
- Turning on the spot: push the handrims in opposite directions at the same time.
- Braking: reduce the speed of the handrims with your hands.

Carry out the first test drives carefully, until you have got used to the wheelchair and its driving characteristics. Wheelchairs only have limited security against tipping over and slipping. Take extra care when braking, starting off or turning on slopes. Extended leaning out of the wheelchair also increases the risk of tipping over. Tipping protection reduces this risk. For easier obstacle negotiation, the attendant can tip the wheelchair (see Chapter 4.15).

7. Technical data

	Pyro Light Optima	Pyro Light Optima XL
Dimensions and weights		
Total length with legrests (mm):	1100	
Total length without legrests (mm):	840	
Total width (mm):	SW + 200	
Total width, folded (mm):	310	
Total width with drum brake (mm):	SW + 200	
Total width with drum brake, folded (mm):	310	
Total height (mm):	860	
Push handle height (mm):	830 - 960	
Lower leg length (mm):	400 - 480	
Kerb weight (kg):	16,8 (sw 430)	18,5 (sw 520)
Weight of heaviest part (kg):	8,8 (sw 430)	10,3 (sw 520)
Max. payload (kg):	125	
Seat width (mm):	370, 400, 430, 460, 490, 520	520, 550, 580
Seat depth (mm):	415, 440*, 465, 490, 515	450
Front seat height without seat cushion (mm):	450, 475, 500	
Rear seat height without seat cushion (mm):	400, 420, 440, 460*, 480, 500	
Backrest angle (°), (Default*)	0 - 11,5 (4*)	
Backrest angle to the seat surface (°):	100	
Centre of gravity (*Standard position)	-28, -8, 12*, 32	
Backrest height without seat cushion (mm):	400, 420, 440, 460	
Armrest height without seat cushion (mm):	190, 205, 220, 235, 250, 265	
Armrest length (mm):	340	
Armrest width (mm):	50	

	Pyro Light Optima	Pyro Light Optima XL
Turning area (mm):	1160	
Static stability (°):	7	
Parking brake operating force (N)	< 60	
Tyres		
Drive wheels (inch):	24" x 1 3/8 PU	
Castors (mm):	8" x 1,25" PU	8" x 2" PU
Material		
Frame:	Aluminium	
Seat and backrest:	Nylon	
Storage conditions		
Temperature:	-25 - +65°C	
Humidity:	45 - 60%	
Operating conditions		
Temperature:	-25 - +50°C	
Humidity:	20 - 80%	
Accessories: Drum brake, anti-tip rollers, lap belt, therapy table, togo, crutch holder, spokeguard, headrest, one-hand operation (wheels), one-hand brake, handrim cover, infusion, oxygen bottle holder.		

The dimensions stated are subject to manufacturing tolerances of +/- 10 mm.

Pyro Light Optima | Lightweight wheelchair

Pyro Light Optima XL

8. Transport

8.1. Transporting the wheelchair

In order to make the wheelchair as manageable as possible to transport, remove all removable parts (legrests, drive wheels, armrests, loose accessories) and fold the wheelchair up (see chapter 4).

8.2. Transportation in vehicles


To see whether your wheelchair is permitted for use as a seat in a Vehicle for Transporting Individuals with Impaired Mobility or not, you can refer to the sticker on the identification label on the wheelchair.

- Wheelchair is designed for use as a seat in a Vehicle for Transporting Individuals with Impaired Mobility (Fig. A).
- Wheelchair is not designed for use as a seat in a Vehicle for Transporting Individuals with Impaired Mobility (Fig. B)

If your wheelchair is not designed to be used as a seat in a vehicle, under no circumstances must it be used as such. It does not fulfil the requirements in compliance with ISO 7176-19. Violation can cause the most serious injuries resulting in death in the case of an accident (Fig. B).

If your wheelchair is not crash-tested, under no circumstances must it be used as a vehicle seat in a vehicle for transporting individuals with impaired mobility. There is a risk of death!

9. Maintenance instructions for the user


9.1. Cleaning

The seat and back materials can be cleaned with warm water using a sponge and/or a soft brush. In case of stubborn dirt, a commercial cleaning agent can be added to the water.

Note: Never use aggressive cleaning agents such as solvents / scouring agents or hard brushes.

- Plastic parts are best cleaned with a commercially available plastic cleaner. Please observe the product information of the cleaning agent.
- Plastic claddings are attacked by non-ionic surfactants and solvents, especially alcohols!
- Wipe the frame parts of the wheelchair with a damp cloth. For heavier soiling, also use a mild cleaning agent.
- Clean the wheels with a damp brush with plastic bristles (do not use a wire brush!).
- Do not use a high-pressure cleaner.
- Wash the push handles, arm pads, seat and back pads with a mild detergent.
 - Upholstery parts can be machine washed at 40°C. Do not machine dry, only dry at room temperature.

9.2. Disinfection

Before disinfection, the wheelchair must be thoroughly cleaned. The wheelchair can be disinfected with a household disinfectant. Use approved spray disinfectants or alcohol-based wipe-on disinfectants.

- In the case of multiple product users, a commercially available disinfectant is prescribed.
- Use only wipe-on disinfectants approved by the Robert Koch Institute (RKI) and the Association for Applied Hygiene (Verbund für Angewandte Hygiene e.V.) (e.g. Bacillol AF or Kohrsolin FF); see manufacturer's recommendations for exposure time and concentration.
- Caution: Damaged upholstery parts cannot be sufficiently disinfected. These should be replaced.
- Under certain circumstances, the surfaces can be affected by disinfectants. This may result in long-term functional limitations. Please note our enclosed overview of the most important care, hygiene and maintenance work.

Attention: Danger of infection! If there are several users, the wheelchair must be cleaned after each use.

Attention! Do not use a high-pressure or steam cleaner.

Please also observe the Bischoff & Bischoff care, warranty, hygiene and maintenance instructions. Your specialist dealer will be happy to provide you with these.

9.3. Checks before driving

9.3.1. Checking the brakes

Before setting off each time, check that the brakes work properly. With the parking brakes applied, the wheelchair cannot be pushed on dry surfaces. When the brakes are released, the wheelchair runs without scraping noises and good directional stability. If this is not the case, notify your dealer immediately and stop using the product until it has been serviced.

9.3.2. Checking the tyres

Each time before setting off, check the tyres for damage.

10. Repairs

Repairs and servicing work must only be carried out by authorised partners (medical supply stores or dealers). To find a partner in your area, please contact our Customer Service. You can find the contact details on the back cover of this user manual.

11. Maintenance information


Typical wear parts are brakes and tyres. Please contact your dealer for spare parts.

We recommend that you have your dealer check your wheelchair at least once a year.

If there are faults or defects, pass the wheelchair to your dealer for servicing immediately.

On request, any information and documents necessary for the repair and servicing will be provided to the dealer. If you have any questions please contact Bischoff & Bischoff Customer Service.

You can find the contact details on the back cover of this user manual.

The wheelchair must be checked in accordance with the maintenance plan and be in proper technical condition.

12. Passing the wheelchair on to another user

The wheelchair is designed for re-use. Before it is passed on to another user, it must be serviced by a dealer and hygienically cleaned. When passing the wheelchair on to another user please remember to hand over this user manual and the address of your dealer to the new user.

13. Storage / shipping

If the wheelchair is to be stored or shipped, all loose parts which are not fixed, except for the drive wheels, must be removed and packed individually in an appropriate cardboard box. The individually-packed parts can then be packed together in a larger cardboard box. We recommend that you keep the original packing and store it for these purposes, so that it is available when needed. Your wheelchair is then

optimally protected against environmental influences during storage or transportation.

For the storage of your wheelchair, please observe the information in the technical data. To avoid the formation of mould and damage to the upholstered parts, the storage location should be as dry as possible and not exposed to direct sunlight.

14. Disposal

If your wheelchair is no longer in use and needs to be disposed of, contact your dealer.

If you want to handle the disposal yourself, find out about the regulations in your area from local recycling companies.

15. Warranty

1. For products supplied by Bischoff & Bischoff, the period of warranty is 24 months from purchase date. The products are delivered free from production and material defects. Where it can be demonstrated that there is a material or manufacturing defect, defective parts will be replaced free of charge.

2. Claims from this warranty shall not arise, if a repair or replacement of a product or a part is required for the following reasons:

- a) normal wear and tear, which include but is not limited to the following parts where fitted: Batteries, motor brushes, handles, armrests, upholstery, tyres, brakes, caps, etc.
- b) overloading the product such as, for example, by exceeding the maximum user weight or load.
- c) the product or part has not been maintained or serviced in accordance with the operating manual, the care and hygiene instructions or the manufacturer's recommendations as shown in the maintenance instructions.
- d) Accessories have been used which are not original accessories.
- e) The product or part has been damaged through negligence, accident or incorrect use.
- f) changes/modifications have been made to the product or parts, which deviate from the manufacturer's specifications.
- g) Repairs have been carried out, before our customer service has been informed of the circumstances.

3. To enforce the warranty, please contact your dealer with the exact details of the nature of the difficulty. The product must be repaired by a dealer authorised by Bischoff & Bischoff.

4. For parts, which have been repaired or exchanged within the scope of the warranty, the warranty will be extended to the remaining warranty period for the product in accordance with point 1).

5. For original spare parts which have been fitted at the customer's expense, these will have a 12 months guarantee, following the fitting.

6. If the rectification of defects fails after a reasonable period of time, the purchaser may, at its discretion, demand a reduction of the remuneration or withdraw from the contract..

7. The warranty is subject to the law of the country in which the Bischoff & Bischoff product was purchased.

In addition, please observe the Bischoff & Bischoff care, warranty, hygiene and maintenance instructions. Your specialist dealer will be happy to provide you with these.

Note!

Despite thorough research and editing, incorrect information may occur in our user manual. We reserve the right to make technical modifications. Images may show accessories. You can find the current version of the User manuals in the download area of our website:

www.bischoff-bischoff.com

Warranty certificate

Product: Pyro Light Optima / Pyro Light Optima XL

Serial number:*

* (To be filled in by the dealer)

Dealer:

Date and stamp

(Pay attention to our general terms and conditions.)

*Signs are located on the cross brace

Failure to observe the instructions in the User manual as well as not carrying out maintenance work properly as well as, in particular, making technical modifications and additions (attachments) without the agreement of Bischoff & Bischoff GmbH will lead to loss of warranty and the manufacturer's liability.

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