
HARTMANN

# Pütter<sup>®</sup>-haft

Kohäsive Kurzzugbinde  
 Cohesive short-stretch bandage  
 Bande cohésive à allongement court  
 Cohesive zwachtel met korte rek  
 Benda coesiva a corta trazione  
 Venda cohesiva de corta elasticidad  
 Ligadura coesiva de curta tração  
 Αυτοσυγκρατούμενος επίδεσμος μικρής έκτασης  
 Kohezyjna opaska o niskiej rozciągliwości  
 Kohezív, rövid megnyúlású pólya  
 Когезионный бинт короткой растяжимости  
 Fașă coezivă cu tracțiune scurtă  
 自黏低拉力繃帶  
 رباط لاصق قصير المدى

1

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## EN Instructions for use

### General description of the device

The short-stretch compression bandages are wide woven. The bandages themselves are made of 100 % cotton and are brown in colour. They have a cohesive layer. The non-active and non-sterile devices are available in different widths and lengths. Pütter-haft<sup>®</sup> are textile-elastic short-stretch bandages.

### Properties

Extensibility: approx. 60 %.

### Intended purpose and indications

The compression bandages are non-sterile medical devices used by healthcare professionals on intact skin of humans (above 6 months of age) in clinical and home environment. The compression bandages can be used for compression therapy (ABPI 0.9-1.3) (decongestion phase) for the treatment of:

- varicose veins (e.g. varicose in pregnancy, supportive in the invasive treatment of varicosis)
- chronic venous insufficiency (CVI) with or without venous leg ulcers (VLU) (e.g. CVI according to CEAP classification C3-C6, primary and secondary prevention of leg ulcer venosum, venous insufficiency and angiodysplasia)
- thrombosis/thrombophlebitis (e.g. thrombophlebitis as well as condition after healed thrombophlebitis, deep leg and arm vein thrombosis, condition after deep vein thrombosis, post-thrombotic syndrome, thromboprophylaxis)
- chronic oedemata of other origin (e.g. oedema in pregnancy, post-traumatic oedema, hormonal oedema, lipoedema, stasis conditions due to immobility (arthrogenic congestion syndrome, paresis), medically induced oedema [e.g. calcium antagonists, isosorbide dinitrate, lithium ointment, sex hormones])
- lymph oedemata

The devices can also be used as:

- support and relief bandages in musculoskeletal injuries

The compression bandages can also be used in combination with primary and secondary dressings. Pütter-haft are single-use devices.

### Composition

100 % cotton

### Contraindications

The short-stretch compression bandages must not be used in case of:

- advanced peripheral arterial occlusive disease (PAOD)
- decompensated heart failure (NYHA III + IV)
- phlegmasia coerulea dolens
- ankle-brachial pressure index (ABPI) > 1.3 and ≤ 0.5

- septic phlebitis
- florid erysipelas

#### Precautions

Special precautionary measures are required in the event of:

- pronounced sensibility disorders of the extremities (e.g. in the case of advanced diabetic peripheral neuropathy)
- pronounced weeping dermatoses
- incompatibility with bandage material
- primary chronic polyarthritis

The application technique, e.g. with regards to stretch and/or padding, may have to be adapted to the patient's anatomical and pathophysiological characteristics, especially in case of fragile or aged skin.

If the foot pulse is not palpable, a Doppler ultrasound should be carried out to determine the ABPI prior to starting therapy. The right measurement method must be considered for the ABPI. During the use of the compression bandages, the ABPI must be controlled regularly. The therapy must be adapted accordingly.

If an inelastic system is created by close to full stretch application, high ankle pressures should be closely monitored by a physician.

The devices should not be used on infants below the age of 6 months. When applying the compression bandages on children (from 6 months on), the pressure should be reduced appropriately. If at any time during therapy the patient develops severe pain, numbness, tingling or noticeable discolouration of the toes, a healthcare professional must be contacted immediately, and the compression bandages must be removed.

#### Application (see pictograms)

Modified bandaging technique according to Pütter, if not otherwise prescribed by the doctor. (1)

#### Further remarks

The devices are indicated for use by healthcare professionals who have received appropriate training. Before application, clean the wound and cover it with a suitable wound dressing. If required, excess padding bandage material may be used for additional padding of bony prominences and tissue recesses (e.g. shinbone, Bisgaard's backdrop). If additional fixation is required, use tape but no bandage clips. After application, the compression bandage must be checked (e.g. circulation in the patient's toes, pain reported by the patient) in order to be certain that the pressure of the bandage is not too high.

Reusing a single-use medical device is dangerous. Re-processing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

#### Incident reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745 on Medical Devices), if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

#### Product disposal

In order to minimize the risk of potential infection hazards, or environmental pollution, disposable components of the medical device should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards. Dispose of the medical device with regular hospital waste.

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**PAUL HARTMANN AG**