MELAcontrol® Helix

Gebrauchsanweisung
Instructions for use
Инструкции за употреба
Brugsanvisning
Kasutusjuhend
Οδηγίες χρήσης
Käyttöohjeilla
Notice d'utilisation
Istruzioni per l'uso
Naudojimo instrukcija
Gebruiksaanwijzing
Instrukcja używania
Instruções de utilização
Instrucțiuni de utilizare
Bruksanvisning
Návod na použitie
Navodila za uporabo
Instrucciones de uso
Návodem k použití
Használati útmutató

MELAcontrol® Helix

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E-Mail: info@melag.de

Web: www.melag.com

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MELAcontrol® Helix

Instructions for use

Email: info@melag.de
Web: www.melag.com

Original instructions

Responsible for content: MELAG Medizintechnik GmbH & Co. KG
We reserve the right to technical alterations

Product name

MELAcontrol Helix - Test body system in accordance with EN 867-5 and EN ISO 11140-1

Intended use

MELAcontrol Helix is a test body system (Process Challenge Device = PCD) for products with narrow lumen. The test body system is used to demonstrate the rapid and complete penetration of steam into medical devices that are sterilised in a steam sterilizer with type "B" cycles. The chemical indicator contained in the sealing cap is a type 2 indicator according to EN ISO 11140-1. If the indicator is exposed to a defined combination of time (min. 3:30 min), steam and temperature (min. 134 °C), it shows a defined colour change from yellow to dark blue/violet. When all air has been removed from the sterilization chamber, steam penetrates the test body system and the indicator shows a uniform colour change. If air or non-condensable gases are present the indicator shows distinct yellow markings.

MELAcontrol Helix supplements the physical measurements and is only suitable for testing the following steam sterilizers:

- steam sterilizers with type "B" cycles in accordance with EN 13060
- steam sterilizers with type "S" cycles in accordance with EN 13060, of which the scope of supply covers the reprocessing of products with narrow lumen
- large steam sterilizers in accordance with EN 285

Hollow body load

Suitable for steam sterilization at 134-137 °C

Label as medical device

CE marking

General guidelines

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

Application

PLEASE NOTE

Only use dry test body systems and, if possible, those that have been blown out with compressed air. Residue condensate in the test body has a negative effect on air removal.

- Insert an unused indicator strip in the plastic bracket of the sealing cap before every test.

The indicator strips are approved for single use only.

1 Unscrew the test body sealing cap.

2 Fold the indicator strip in half.

3 Insert the indicator strip, fold first, as far as possible into the sealing cap.

4 Screw the sealing cap hand-tight onto the test body.

**CAUTION**

Danger of burns from hot condensate!

Position the test body for sterilization in a level position, otherwise hot condensate can run out of the test body after removal.

- Let the test body cool before removal.

Use for daily routine check

1. Place a new indicator strip in the test body.
2. Slide the test body with the flat side facing downwards onto a tray in the empty sterilization chamber.
3. Start the steam sterilizer in a validated sterilization program according to the intended use.
4. Remove the test body from the steam sterilizer when the program run has ended.
5. Remove the indicator strip from the sealing cap. Evaluate the indicator strip.

Use as batch control

1. Place a new indicator strip in the test body.
2. Place the test body in the sterilization chamber above the load, with the flat side facing down.
3. Start the steam sterilizer in a validated sterilization program according to the intended use.
4. Remove the test body from the steam sterilizer when the program run has ended.
5. Remove the indicator strip from the sealing cap. Evaluate the indicator strip.

Evaluation

	completely changed indicator strip	Sterilisation successful sufficient air removal and steam penetration
	changed indicator strip with slight gradient to the middle	Sterilisation not successful only temperature impact, no steam penetration
	incompletely changed indicator strip	Sterilisation not successful insufficient holding time and steam penetration

Documentation

The evaluation of the results must be performed by a sufficiently trained person holding responsibility for this task. The colour change is caused by an irreversible chemical reaction. If stored correctly, the indicator strips can be archived for years. EN ISO 15882 only requires a description of the results. This can be performed in written or digital form. Comply with the valid national and regional regulations.

Important information regarding routine operation

Replace the test body after 250 indicator strips have been used.

Please comply with the current recommendations issued by the Robert Koch Institute (RKI) and the information contained in DIN 58946-7. Further information regarding routine operation is to be found in user manual of the steam sterilizer.

Storage

Store the indicator strips in a cool, dark and dry place before and after use.

Do not store the indicator strips together with chemical substances and gases (e.g. ethylene oxide [EO]), hydrogen peroxide [H₂O₂] or instruments sterilised with H₂O₂. Incorrect storage of the indicator strips can cause unwanted colour changes. Further information regarding storage is provided on the packaging. A secure determination of malfunctions in the process can only be guaranteed up to the specified date of expiry.

Disposal

Dispose of this product and spare parts, e.g. seals, that are no longer used properly. Comply with all relevant disposal specification in terms of possibly contaminated waste.

The packaging protects the product against transport damage. The packaging materials have been selected for their environmentally-friendly and recycling properties and can be recycled. Returning the packaging to the material flow reduces the amount of waste and saves raw materials.

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