

Product Name	i-gel® Plus, supraglottic airway		
Reference numbers	8603000	8604000	8605000
Product Type	i-gel® supraglottic airway		
Manufacturer	Intersurgical		
Materials		ABS	PP
		TPE	
Packaging			
Unit Packaging	Individually packaged		
Box quantity	10		
Pallet quantity	640		
Instructions for Use			
Instructions for use with product	No		
Instructions for use with each box	Yes		
CE Class	IIa		
Applicable Standards	BS EN ISO 13485: Medical devices — Quality management systems — Requirements for regulatory purposes ISO 9001: Quality Management Systems - Requirements ISO 15223-1: Medical Devices - Symbols to be Used with Medical Device Labels, Labelling, & Information to be Supplied BS EN ISO 5356-1: Anaesthetic and respiratory equipment - Conical connectors. Cones and sockets BS EN 13544-2: Respiratory therapy equipment. Tubing and connectors BS ISO 11712: Anaesthetic and respiratory equipment. Supralaryngeal airways and connectors		
Sterile	Yes		
Shelf life	3 years		
Disposal Methods	As hospital guidelines		
Re-use	 Single use only		

Product characteristics

Technical information

	8603000 Small adult	8604000 Medium adult	8605000 Large adult
i-gel size	3	4	5
Patient weight (kg)	30–60	50–90	90+
Internal volume of airway channel (ml)	20.3	25.4	30.1
Weight including packaging (g)	109	127	150



**Material
Abbreviations**

ABS	Acrylonitrile Butadiene Styrene
AL	Aluminium
CR	Polychloroprene Rubber
EDPM	Ethylene Propylene Terpolymer Rubber
HDPE	High Density Polyethylene
HIPS	High Impact Polystyrene
LDPE	Low Density Polyethylene
LSR	Liquid Silicone Rubber
MBS	Methacrylate Butadiene Styrene
NRL	Natural Rubber Latex
PA	Polyamide (Nylon)
PC	Polycarbonate
PE	Polyethylene
POM	Polyoxymethylene
PP	Polypropylene
PS	Polystyrene
PSU	Polysulfone
PU	Polyurethane
PVC	Poly Vinyl Chloride
SB	Styrene Butadiene
Si	Silicone
SS	Stainless Steel
TPE	Thermoplastic Elastomer
SBC	Styrene Butadiene Copolymer
SEBS	Styrene-Ethylene-Butadiene-Styrene

