


Product Name	i-view video laryngoscope	
Reference numbers	8008000	
Product Type	Video Laryngoscope	
Manufacturer	Intersurgical	
Materials	 Lens cover seal - TPE Lens cover - Acrylic Battery spring - Nickel plated steel Casing - Glass filled nylon Bezel & Button - ABS	
Packaging		
Unit Packaging	Individually packaged	
Box quantity	8	
Pallet quantity	512	
Instructions for Use	Instructions for use with product Instructions for use with each box NO YES	
CE Class	I	
Applicable Standards	BS EN ISO 13485: Medical devices — Quality management systems — Requirements for regulatory purposes ISO 9001: Quality Management Systems - Requirements ISO 15223-1: Medical Devices - Symbols to be Used with Medical Device Labels, Labelling, & Information to be Supplied IEC 60601-1: Medical electrical equipment. General requirements for basic safety and essential performance. IEC 60601-1-6:2010+A1:2013 - Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability IEC 60601-1-12: Medical electrical equipment - Part 1-12: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems intended for use in the emergency medical services environment	
Sterile	No	
Shelf life	5 Years	
Disposal Methods	As hospital guidelines	
Re-use	 Single use only	

Product characteristics

Technical information -

Size (mm)	A - Height	202
	B - Depth	161
	C - Width	70
Weight (g)		141
Electrical rating		4.5v DC 150mA



**Material
Abbreviations**

ABS	Acrylonitrile Butadiene Styrene
AL	Aluminium
CR	Polychloroprene Rubber
EDPM	Ethylene Propylene Terpolymer Rubber
EVA	Ethyl Vinyl Acetate
HDPE	High Density Polyethylene
HIPS	High Impact Polystyrene
LDPE	Low Density Polyethylene
LSR	Liquid Silicone Rubber
MBS	Methacrylate Butadiene Styrene
NRL	Natural Rubber Latex
PA	Polyamide (Nylon)
PC	Polycarbonate
PE	Polyethylene
POM	Polyoxymethylene
PP	Polypropylene
PS	Polystyrene
PSU	Polysulfone
PU	Polyurethane
PVC	Poly Vinyl Chloride
SB	Styrene Butadiene
Si	Silicone
SS	Stainless Steel
TPE	Thermoplastic Elastomer
SBC	Styrene Butadiene Copolymer
SEBS	Styrene-Ethylene-Butadiene-Styrene

