

«BASIC» VAGINAL SPECULA



DESCRIPTION

The disposable **Teqler «Basic» vaginal specula** are used for gynaecological examinations and made from transparent plastic. The single use Cusco specula are available in 3 colour-coded sizes.

PRODUCT DETAILS

- Cusco vaginal specula
- Transparent plastic
- With lock (basic model)
- Available in S, M and L, colour-coded
- No risk of infectious transmissions
- Suitable for examining the vaginal lining
- Individually sterile packed
- 1 pack of 100 specula

SPECULA SIZES

- Size S: 11 x 4 cm (colour: blue)
- Size M: 12,5 x 4,5 cm (colour: red)
- Size L: 13 x 5 cm (colour: green)

Overview

Article Number:
T135390

Category:
**Gynaecological
Instruments**

Teqler is a trademark of NetMed S.à.r.l.

 NetMed S.à.r.l.
8, Am Scheerleck, 6868 Wecker, Luxembourg

 +352 267149 09  +352 267149 19  info@netmed.lu

Social Media

