

Instruction for use

Cardboard Mouthpiece

Shanghai Everpure Medical Plastic Co., Ltd.

1. Indications

Cardboard Mouthpiece is used to fix patient's mouth during the spirometer testing..

2. Product specification(s) and description

Model:

Description: Cardboard Mouthpiece consists of cardboard.

3. Instruction for use

Spirometer Testing:

- 1) Take off the previous cardboard mouthpiece
- 2) Insert a new cardboard mouthpiece into the spirometer
- 3) Let patient keep the mouthpiece in the mouth and blow the air into the spirometer.

4. Precautions

Please do replace the cardboard mouthpiece every time.

5. Sterilization

EO Sterilization

6. Shelf Life

5 years

7. Label

- Symbol  for "DO NOT REUSE" (The product is for single use.)
- Symbol  for "BATCH CODE" (the symbol should be attached with batch code and be close to the graph, batch code, lot number and batch number). Examples: LOT ABC123
- Symbol  for "USE BY" like YYYYMMDD
- Symbol  for "DATE OF MANUFACTURE" (the symbol should be close to the number)
- Symbol  for "MANUFACTURE"

- Symbol  for "EC REPRESENTATIVE"
- Symbol  for "ATTENTION, SEE INSTRUCTIONS FOR USE"
- Symbol  for "Symbol for "CONSULT INSTRUCTIONS FOR USE"
- Symbol  for "TEMPERATURE LIMITATION"
- Symbol  for " Keep away from sunlight "
- Symbol  for " Avoid rain "
- Symbol  for " CE " Marking
- Symbol  for " CONTAINS OR PRESENCE OF PHTHALATE"
- Symbol  for " Thrown into the trash after use "
- Symbol  for "Do not use if package is damaged"
- Symbol  for "The product is sterilized by EO"

Manufacturer ☐

Shanghai Everpure Medical Plastic Co., Ltd.

Address: No.2418 Tingfeng Rd., Xinnong Town,
Jinshan District, Shanghai 201503, China

Tel/Fax: 0086-21-57340011, 57340022